

# TESTAmed COVID-19 Antigen Rapid Test

## Lateral Flow Chromatographic Immunoassay

### Performance Summary



# Table of Contents

<b>TESTAmed COVID-19 Antigen Rapid Test</b> .....	1
Introduction .....	1
Key Features.....	1
Objective .....	1
Product Specifications.....	2
Instruction for Use .....	3
Interpretation of Results .....	3
<b>Performance Evaluation</b> .....	4
Clinical Sensitivity and Specificity Evaluation .....	4
Analytical Sensitivity – Limit of Detection (LoD).....	6
Analytical Specificity – Cross-Reactivity Study .....	9
Interference Studies .....	11
<b>Summary</b> .....	13

# TESTAmed COVID-19 Antigen Rapid Test

## Introduction

The coronavirus disease 2019 (COVID-19) is caused by the severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2). With the spread of the COVID-19 pandemic, the demand for tools to quickly scan for the disease has rapidly increased worldwide.

The TESTAmed COVID-19 Antigen Rapid Test is a lateral flow immunoassay for the qualitative detection of nucleocapsid protein antigen from SARS-CoV-2. The sample material is obtained by using a nasopharyngeal swab during the acute phase of the infection. The virus detection takes place in 15 minutes, for quick and clear results.

The antigen test is confirmed to be sensitive to the new COVID-19 variants with spike protein mutations discovered in the UK and South Africa.

## Key Features

- Accurate and Reliable – clinical study shows 94.2% sensitivity and 99.6% specificity
- Convenient – complete a test in 5 simple steps
- Rapid – results in only 15 minutes

## Objective

This document was created to evaluate the clinical performance and reliability of TESTAmed COVID-19 Antigen Rapid Test according to clinical and laboratory standards. The goal is to demonstrate the excellent performance of this Lateral Flow Chromatographic Immunoassay.








## Product Specifications

Table 1 – TESTAmed COVID-19 Antigen Rapid Test Specifications





Test Principle	Lateral Flow Chromatographic Immunoassay
Target Antigen	SARS-CoV-2 Nucleocapsid Protein
Sample Type	Fresh Nasopharyngeal Swab Specimen
Limit of Detection (LoD)	1.26 x 10 <sup>2</sup> TCID <sub>50</sub> per mL
Cross-Reactivity & Interferences	Viruses, Bacteria and Interferences tested do not cross-react or interfere
Reaction time	15 minutes



## Instruction for Use

<p><b>1</b></p>  <p>Add 10 drops (500 <math>\mu</math>L) of extraction buffer into the extraction tube.</p>	<p><b>2</b></p>  <p>10~20 times</p> <p>Immerse the patient nasopharyngeal swab sample into the extraction tube. Roll the swab at least 10 - 20 times while pressing the head against the bottom and side of the Extraction Tube.</p>	<p><b>3</b></p>  <p>30 sec</p> <p>Leave the swab in the extraction buffer for 30 seconds.</p>	<p><b>4</b></p>  <p>5~10 times</p> <p>When removing, roll the swab head toward the inside of the extraction tube and squeeze the sides of the tube to extract the liquid from the swab. Dispose of the used swab in your biohazardous waste.</p>
<p><b>5</b></p>  <p>Place the nozzle cap and press tightly.</p>	<p><b>6</b></p>  <p>3~5 drops</p> <p>Add 3-5 drops (100 <math>\mu</math>L) of the processed sample into the sample well.</p>	<p><b>7</b></p>  <p>15-20 min</p> <p>Visual interpretation within 15-20 minutes. Some positive results may appear sooner.</p> <p><b>CAUTION:</b> Do not read the results after 20 minutes. It may provide false results.</p>	

## Interpretation of Results

<p>Positive</p>		<p><b>+</b> <b>Positive result</b></p> <p>There are two colored lines on the test cassette. Both colored test and control lines appear on the test cassette. Within the specified observation time, a weak colored test line should be judged as a positive result.</p>
<p>Negative</p>		<p><b>-</b> <b>Negative result</b></p> <p>Only the colored control line appears on the test cassette. The absence of the test line indicates a negative result.</p>
<p>Invalid Assay</p>		<p><b>×</b> <b>Invalid result</b></p> <p>There should always be a colored control line in the control region regardless of the test result. If the control line is not seen, repeat the assay with a new test cassette.</p>
<p>Negative / High CT value</p>		<p><b>Negative / High CT value</b></p> <p>Within a specified observation time, a very weak color on the T line should be judged as negative or high CT value. Please refer to the RT-PCR result.</p>

# Performance Evaluation

## Clinical Sensitivity and Specificity Evaluation

### Objective

This clinical evaluation study is intended to evaluate the clinical sensitivity and clinical specificity of COVID-19 Antigen Rapid Test by comparing the test results to an RT-PCR test.

Clinical sensitivity and specificity of fresh nasopharyngeal swab specimens were measured by healthcare professionals or non-laboratory personnel in point-of-care settings or near-patient sites.

### Method

Clinical sensitivity and specificity of fresh nasopharyngeal swab specimens were measured by healthcare professionals or non-laboratory personnel in point-of-care settings or near-patient sites. The study collected samples from subjects from Taiwan, USA and Bosnia with positive (+) and negative (-) results.

Table 2 – Study Site Subjects

Site	COVID-19 (+)	COVID-19 (-)	Total
Taiwan	36 subjects	100 subjects	136 subjects
USA	7 subjects	138 subjects	145 subjects
Bosnia	60 subjects	30 subjects	90 subjects
Total	103 subjects	268 subjects	371 subjects

### Acceptance Criteria

Sensitivity = Positive Percent Agreement, PPA (%)  $\geq$  80% and  $\geq$  90% (Ct<30)

Specificity = Negative Percent Agreement, NPA (%)  $\geq$  97%

## Results

### Overall

Test Device: TESTAmed COVID-19 Antigen Rapid Test

Comparator Method: RT-PCR test

Table 3 – Clinical Performance of TESTAmed COVID-19 Antigen Rapid Test

		PCR Test Results		
		POS	NEG	Subtotal
TESTAmed COVID-19 Antigen Rapid Test	Positive	97	1	98
	Negative	6	267	273
	Subtotal	103	268	371
Sensitivity		94.2% (95%CI: 87.9%-97.3%)		
Specificity		99.6% (95%CI: 97.9%-99.9%)		

## Conclusion

The clinical performance of TESTAmed COVID-19 Antigen Rapid Test was determined by testing 103 positive and 268 negative specimens for SARS CoV-2 antigen (Ag). Based on the clinical study, the antigen test has a sensitivity of 94.2% (95%CI: 87.9%-97.3%) and specificity 99.6% (95%CI: 97.9%-99.9%)The results meet the acceptance criteria of sensitivity  $\geq 80\%$  and  $\geq 90\%$  (Ct<30), and specificity  $\geq 97\%$ .

## **Analytical Sensitivity – Limit of Detection (LoD)**

### **Objective**

This study aims to determine the limit of detection (LoD) of the TESTAmed COVID-19 Antigen Rapid Test, and to simulate the detection method when using nasopharyngeal (NP) swabs.

### **Method**

The Limit of Detection (LoD) of TESTAmed COVID-19 Antigen Rapid Test was determined using limiting dilutions of live SARS-CoV-2, isolate TWN/CGMH-CGU-01. In this study, all the SARS-CoV-2 serial dilutions were made in the SARS-CoV-2 negative nasopharyngeal swab pool. The material was supplied frozen at a concentration of  $10^{5.4}$  TCID<sub>50</sub> per mL. The study to determine the TESTAmed COVID-19 Antigen Rapid Test LoD was designed to reflect the assay when using direct nasopharyngeal swab.

The LoD was determined in three steps:

- 1) LoD Screening: 10-fold dilutions of the live SARS-CoV-2 were made. These dilutions were tested in triplicate. The concentration demonstrating 3 of 3 positives was chosen for LoD range finding.
- 2) LoD Range Finding: Five (5) 2-fold dilutions of the suitable concentration were made. These dilutions were tested in triplicate. The concentration demonstrating 3 of 3 positives was chosen for LoD confirmation.
- 3) LoD Confirmation: The dilution was tested for a total of twenty (20) results for LOD confirmation.

### **Acceptance Criteria**

The LoD concentration must be  $< 10^3$  TCID<sub>50</sub>/mL.



## Results

Table 4 – The results of LoD Screening

TESTAmed COVID-19 Antigen Rapid Test	Negative human clinical NP Swab pooled extraction	Spiked virus concentration (TCID <sub>50</sub> /ml)					
	0	1	1/10	1/100	1/1000	1/10000	1/10000
	0	10 <sup>5.4</sup>	10 <sup>4.4</sup>	10 <sup>3.4</sup>	10 <sup>2.4</sup>	10 <sup>1.4</sup>	2.5
Test 1	-	+	+	+	+	-	-
Test 2	-	+	+	+	+	-	-
Test 3	-	+	+	+	+	-	-
Overall call rate	0/3 (0%)	3/3 (100%)	3/3 (100%)	3/3 (100%)	3/3 (100%)	0/3 (0%)	0/3 (0%)

Table 5 – The results of LoD Range Finding

TESTAmed COVID-19 Antigen Rapid Test	Negative human clinical NP Swab pooled extraction	Spiked virus concentration (TCID <sub>50</sub> /ml)					
	0	1/1000	1/2000	1/4000	1/8000	1/16000	1/32000
	0	10 <sup>2.4</sup>	1.26x10 <sup>2</sup>	0.63x10 <sup>2</sup>	0.31x10 <sup>2</sup>	0.16x10 <sup>2</sup>	0.08x10 <sup>2</sup>
Test 1	-	+	+	-	-	-	-
Test 2	-	+	+	-	-	-	-
Test 3	-	+	+	-	-	-	-
Overall call rate	0/3 (0%)	3/3 (100%)	3/3 (100%)	0/3 (0%)	0/3 (0%)	0/3 (0%)	0/3 (0%)

Table 6 – The results of LoD confirmation

TESTAmed COVID-19 Antigen Rapid Test	Negative human clinical NP Swab pooled extraction	Spiked virus concentration (TCID <sub>50</sub> /ml)		
		1/1000	1/2000	1/4000
	0	10 <sup>2.4</sup>	1.26x10 <sup>2</sup>	0.63x10 <sup>2</sup>
Test 1	-	+	+	-
Test 2	-	+	+	-
Test 3	-	+	+	-
Test 4	-	+	+	-
Test 5	-	+	+	-
Test 6	-	+	+	-
Test 7	-	+	+	-
Test 8	-	+	+	-
Test 9	-	+	+	-
Test 10	-	+	+	-
Test 11	-	+	+	-
Test 12	-	+	+	-
Test 13	-	+	+	-
Test 14	-	+	+	-
Test 15	-	+	+	-
Test 16	-	+	+	-
Test 17	-	+	+	-
Test 18	-	+	+	-
Test 19	-	+	+	-
Test 20	-	+	+	-
Overall call rate	0/20 (0%)	20/20 (100%)	20/20 (100%)	0/20 (0%)

## Conclusion

Based on this test, the TESTAmed COVID-19 Antigen Rapid Test concentration of LoD was confirmed as 1.26 x 10<sup>2</sup> TCID<sub>50</sub> per mL

# Analytical Specificity – Cross-Reactivity Study

## Objective

This study aims to demonstrate that the TESTAmed COVID-19 Antigen Rapid Test will not cross-react or interfere with the medically relevant levels of viruses and bacteria listed in the *results* section.

## Method

Cross-reactivity of the TESTAmed COVID-19 Antigen Rapid Test was evaluated by testing various viruses (16) and bacteria (18). Each virus or bacteria was tested in triplicate in the absence or presence of  $3.78 \times 10^2$  TCID<sub>50</sub>/mL (3 LoD) of live SARS-CoV-2. The final concentration of each virus or bacteria was listed in the Table below. Testing was performed in triplicate.

## Acceptance Criteria

The tested concentration of viruses or bacteria should not cross-react or interfere with the results of TESTAmed COVID-19 Antigen Rapid Test.

## Results

Table 7 – Cross Reactivity and Interference of Virus and Bacteria

Virus/Bacteria	Concentration	Cross-Reactive Results	SARS-CoV-2 Concentration (3 LoD)	Interference Results
Human Coronavirus OC43	$2.5 \times 10^5$ pfu/mL	Negative	$3.78 \times 10^2$ TCID <sub>50</sub> /mL	Positive
Human Coronavirus 229E	$2.5 \times 10^5$ pfu/mL	Negative	$3.78 \times 10^2$ TCID <sub>50</sub> /mL	Positive
Influenza A · H1N1	$2.5 \times 10^5$ pfu/mL	Negative	$3.78 \times 10^2$ TCID <sub>50</sub> /mL	Positive
Influenza A · H3N2	$2.5 \times 10^5$ pfu/mL	Negative	$3.78 \times 10^2$ TCID <sub>50</sub> /mL	Positive
Influenza B · Victoria	$2.5 \times 10^5$ pfu/mL	Negative	$3.78 \times 10^2$ TCID <sub>50</sub> /mL	Positive
Influenza B · Yamagata	$2.5 \times 10^5$ pfu/mL	Negative	$3.78 \times 10^2$ TCID <sub>50</sub> /mL	Positive
Respiratory syncytial virus	$2.5 \times 10^5$ pfu/mL	Negative	$3.78 \times 10^2$ TCID <sub>50</sub> /mL	Positive
Rhinovirus	$2.5 \times 10^5$ pfu/mL	Negative	$3.78 \times 10^2$ TCID <sub>50</sub> /mL	Positive
Adenovirus type 1 (Adenoid 71)	$2.5 \times 10^5$ pfu/mL	Negative	$3.78 \times 10^2$ TCID <sub>50</sub> /mL	Positive
Adenovirus type 7	$2.5 \times 10^5$ pfu/mL	Negative	$3.78 \times 10^2$ TCID <sub>50</sub> /mL	Positive
Enterovirus 68	$2.5 \times 10^5$ pfu/mL	Negative	$3.78 \times 10^2$ TCID <sub>50</sub> /mL	Positive
Human parainfluenza type 1	$2.5 \times 10^5$ pfu/mL	Negative	$3.78 \times 10^2$ TCID <sub>50</sub> /mL	Positive

Human parainfluenza type 2	2.5 x 10 <sup>5</sup> pfu/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
Human parainfluenza type 3	2.5 x 10 <sup>5</sup> pfu/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
Human parainfluenza type 4	2.5 x 10 <sup>5</sup> pfu/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
Respiratory syncytial virus type B	2.5 x 10 <sup>5</sup> pfu/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
<i>Bordetella pertussis</i>	2 x 10 <sup>6</sup> CFU/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
<i>Chlamydia pneumoniae</i>	2 x 10 <sup>6</sup> CFU/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
<i>Corynebacterium sp.</i>	2 x 10 <sup>6</sup> CFU/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
<i>Escherichia coli</i>	2 x 10 <sup>6</sup> CFU/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
<i>Hemophilus influenzae</i>	2 x 10 <sup>6</sup> CFU/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
<i>Lactobacillus sp.</i>	2 x 10 <sup>6</sup> CFU/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
<i>Moraxella catarrhalis</i>	2 x 10 <sup>6</sup> CFU/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
<i>Mycobacterium tuberculosis (avirulent)</i>	2 x 10 <sup>6</sup> CFU/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
<i>Neisseria meningitidis</i>	2 x 10 <sup>6</sup> CFU/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
<i>Neisseria sp.</i>	2 x 10 <sup>6</sup> CFU/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
<i>Pseudomonas aeruginosa</i>	2 x 10 <sup>6</sup> CFU/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
<i>Staphylococcus aureus (Protein A producer)</i>	2 x 10 <sup>6</sup> CFU/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
<i>Staphylococcus epidermidis</i>	2 x 10 <sup>6</sup> CFU/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
<i>Streptococcus pneumoniae</i>	2 x 10 <sup>6</sup> CFU/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
<i>Streptococcus pyogenes</i>	2 x 10 <sup>6</sup> CFU/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
<i>Streptococcus salivarius</i>	2 x 10 <sup>6</sup> CFU/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
Pooled human nasal wash – representative of normal respiratory microbial flora	2 x 10 <sup>6</sup> CFU/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
<i>Mycoplasma pneumoniae</i>	2 x 10 <sup>6</sup> CFU/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive

## Conclusion

Cross-reactivity of the TESTAmed COVID-19 Antigen Rapid Test was evaluated by testing various viruses (16) and bacteria (18). Each virus or bacteria was tested in triplicate in the absence or presence of 3.78 x 10<sup>2</sup> TCID<sub>50</sub>/mL (3 LoD) of live SARS-CoV-2.

Based on the data generated by this study, the tested viruses and bacteria do not cross-react or interfere with TESTAmed COVID-19 Antigen Rapid Test results.

## Interference Studies

### Objective

This study aims to demonstrate that the TESTAmed COVID-19 Antigen Rapid Test results will not be interfered with 20 potential interfering substances and concentrations listed in the *results* section.

### Method

Twenty (20) potentially interfering substances that may be found in the upper respiratory tract were selected and tested against the TESTAmed COVID-19 Antigen Rapid Test. Each substance was tested in triplicate in the absence or presence of  $3.78 \times 10^2$  TCID<sub>50</sub>/mL (3 LoD) of live SARS-CoV-2.

### Acceptance Criteria

The tested concentration of interfering substance should not cross-react or interfere with the results of TESTAmed COVID-19 Antigen Rapid Test.

### Results

Table 8 – Cross Reactivity and Interference of Interfering Substances

Interfering Substance	Active Ingredient	Concentration	Cross-Reactive Results	SARS-CoV-2 Concentration (3 LoD)	Interference Results
Ephrine Nasal Spray "GCPC"	Oxymetazoline	5% v/v	Negative	$3.78 \times 10^2$ TCID <sub>50</sub> /mL	Positive
Chloraseptic, Regular strength	Benzocaine / Menthol	1.5 mg/mL	Negative	$3.78 \times 10^2$ TCID <sub>50</sub> /mL	Positive
Tamiflu	Osetamivir	2.5 mg/mL	Negative	$3.78 \times 10^2$ TCID <sub>50</sub> /mL	Positive
Physiomer Saline nasal spray	Saline	15% v/v	Negative	$3.78 \times 10^2$ TCID <sub>50</sub> /mL	Positive
Tobrex Eye Ointment	Tobramycin	51.4 µmol/L	Negative	$3.78 \times 10^2$ TCID <sub>50</sub> /mL	Positive
Sucrets	Dyclonine / Menthol	1.5 mg/mL	Negative	$3.78 \times 10^2$ TCID <sub>50</sub> /mL	Positive
NeilMed NasoGEL Spray	sodium hyaluronate / Saline	5% v/v	Negative	$3.78 \times 10^2$ TCID <sub>50</sub> /mL	Positive
Acetaminophen	Acetaminophen	1324 µmol/L	Negative	$3.78 \times 10^2$ TCID <sub>50</sub> /mL	Positive

Acetylsalicylic acid	Acetylsalicylic acid	3.62 mmol/L	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
Ibuprofen	Ibuprofen	2.425 mmol/L	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
Erythromycin	Erythromycin	81.6 µmol/L	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
Fisherman's Friend	Menthol	1.5 mg/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
Plaquenil	Hydroxychloroquine sulphate	150 µmol/L	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
SUPEROCIN	Ciprofloxacin	30.2 µmol/L	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
Zeffix	Lamivudine	1 mg/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
Blood (human)	Blood (human)	2.5% v/v	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
Ricola	Menthol	1.5 mg/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
Mupirocin	Mupirocin	10 mg/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
Flonase	Fluticasone	5% v/v	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive
Purified mucin protein	Mucin protein	2.5 mg/mL	Negative	3.78 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	Positive

## Conclusion

The interference study of TESTAmed COVID-19 Antigen Rapid Test was evaluated by testing various interfering substances. Each interfering substance was tested in triplicate in the absence or presence of 3.78 x 10<sup>2</sup> TCID<sub>50</sub>/mL (3 LoD) of live SARS-CoV-2.

Based on the data generated by this study, the tested interfering substances do not-cross-react or interfere with TESTAmed COVID-19 Antigen Rapid Test results.

## Summary

The results displayed throughout this report show that the TESTAmed COVID-19 Antigen Rapid Test provides excellent clinical performance results, with sensitivity of 94.2% (positive percent agreement) and specificity of 99.6% (negative percent agreement).

Additionally, the cross-reactivity study and interference study demonstrate the virus, bacteria, and substances tested throughout the study do not interfere with the TESTAmed COVID-19 Antigen Rapid Test results.

## **Performance Summary**

**created by**

 TaiDoc Technology Corp.

B1-7F, No. 127, Wugong 2nd Rd.

Wugu Dist., 24888 New Taipei City

**for**

Distributor:

Sebapharma GmbH & Co. KG

Postfach 1120

56136 Boppard – Bad Salzig

Importer:

Trend Pharma GmbH

Binger Str. 80

56154 Boppard – Bad Salzig

**[www.testamed.de](http://www.testamed.de)**